

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(Benzotriazol-1-yloxy)tripyrrolidinophosphonium hexafluorophosphate - purum, ≥97.0% (TLC)

Product Number: **12805**
CAS Number: 128625-52-5

Formula: C₁₈H₂₈F₆N₆OP₂
Formula Weight: 520.39 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Colour)	White to Light Brown
Appearance (Form)	Powder or Crystals
Purity (HPLC)	≥ 97.0 %
Purity (TLC)	97 - 100 %
Melting Temperature	140 - 160 deg
Proton NMR spectrum	Conforms to Structure

Specification: PRE.0.ZQ5.10000018220

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.