

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
CARBON NANOTUBE ARRAY, MULTI-WALLED, DR&

Product Number: 901046
CAS Number: 308068-56-6

Formula: C

TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Solid
ICP Major Analysis Confirms Silicon Component	Confirmed
Height	0.5 - 1.5 mm
Diameter 10-40 nm	Confirmed
Carbon Content	≥ 98 %

Specification: PRD.0.ZQ5.10000086807

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.