

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,8-Dihydroxynaphthalene - 95%


Product Number: **740683**

CAS Number: 569-42-6

Formula: C₁₀H₈O₂

Formula Weight: 160.17 g/mol

Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

White to Beige and Faint Brown to Brown

Appearance (Form)

Powder or Crystals

Proton NMR Spectrum

Conforms to Structure

Purity (HPLC)

≥ 94.5 %

Specification: PRD.1.ZQ5.10000035742

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.