

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Carbon nanofibers – graphitized (iron-free), composed of conical platelets, D × L 100 nm × 20–200 µm

Product Number: 719781

Formula: C

Formula Weight: 12.01 g/mol

TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Iron (Fe)	≤ 100 ppm

Specification: PRD.0.ZQ5.10000050060

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.