

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

tert-Butyl 2-(4-[[4-(3-azidopropoxy)phenyl]azo]benzamido)ethylcarbamate - 95%

Product Number:

771066

CAS Number:

1207289-71-1

Formula:


C₂₃H₂₉N₇O₄

Formula Weight:

467.52 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

White to Very Dark Yellow and Faint

Orange to Very Dark Orange and Faint Brown to Brown

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Carbon

55.5 - 62.6 %

Nitrogen

19.7 - 22.2 %

Purity (HPLC)

≥ 94.0 %

Specification: PRD.0.ZQ5.10000035526

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.