

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Poly(ethylene glycol) bisazide - average M_n 1,100

Product Number:

756598

CAS Number:

82055-94-5



Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

Colorless and White to Beige and Faint

Brown to Brown

Appearance (Form)

Liquid or Semi-Solid or Paste or Solid

Proton NMR spectrum

Conforms to Structure

Molecular Number

900 - 1300

by HNMR

Poly Dispersity MW/MN

≤ 1.1

Carbon

Nitrogen

Specification: PRD.1.ZQ5.10000020068

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.