


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Furancarboxaldehyde - $\geq 97\%$

Product Number: 278866
CAS Number: 498-60-2
MDL: MFCD00010424
Formula: C₅H₄O₂
Formula Weight: 96.08 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Yellow to Brown
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	$\geq 97.0 \%$

Specification: PRD.1.ZQ5.10000027573

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.