

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Terrific Broth, modified - powder

Product Number: T0918

TEST	Specification
Appearance (Color)	Faint Yellow to Dark Yellow to Dark Beige
Appearance (Form)	Powder
Solubility (Color)	Yellow to Dark Yellow and Brown-Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
1 g plus 10 ml H ₂ O	
pH	7.1 - 7.4
(1% in H ₂ O)	
Suitability	Conforms
for use as microbial culture reagent	

Specification: PRD.1.ZQ5.10000043813

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.