


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-(Dodecylthiocarbonothioylthio)-2-methylpropionic acid pentafluorophenyl ester

Product Number: **740810**
CAS Number: 1174764-26-1



Formula: C₂₃H₃₁F₅O₂S₃
Formula Weight: 530.68 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Very Dark Yellow and Colorless to Very Dark Brown
Appearance (Form)	Liquid or Solid
Proton NMR spectrum	Conforms to Structure
Fluorine NMR Spectra	Conforms to Structure
Purity (HPLC)	≥ 88.0 %

Specification: PRD.2.ZQ5.10000025511

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.