


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Nitrothioanisole - 96%

Product Number: **557471**
CAS Number: 701-57-5
MDL: MFCD00010868
Formula: C₇H₇NO₂S
Formula Weight: 169.20 g/mol



TEST	Specification
Appearance (Color)	Yellow to Yellow-Green
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	47.4 - 51.9 %
Nitrogen	7.9 - 8.7 %
Sulfur	18.1 - 19.8 %

Specification: PRD.0.ZQ5.10000056219

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.