

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name: Anti-Human IgG (whole molecule)-FITC antibody produced in rabbit - IgG fraction of antiserum, buffered aqueous solution  
Product Number: F4512

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	10.0 - 20.0 mg/ml
F/P Molar Ratio	2.5 - 6.5
IF Cells	Hep-2 Human Cells
IF-Titer	1:64-1:128
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000013983

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.