


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Phenyl-5-benzimidazolesulfonic acid - 96%

Product Number: 437166
CAS Number: 27503-81-7
MDL: MFCD00150559
Formula: C₁₃H₁₀N₂O₃S
Formula Weight: 274.3 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	95.5 - 104.5 %

Specification: PRD.0.ZQ5.10000037034

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.