SIGMA-ALDRICH®


sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: 5-Methyl-2-hepten-4-one, predominantly trans - 98%, FG

Product Number: CAS Number: MDL: Formula: Formula Weight: W376108 81925-81-7 MFCD00792509 C8H14O 126.2 g/mol



TEST

Appearance (Color) Appearance (Form) Refractive index at 20 ° C Infrared spectrum Purity (GC) Color Test Arsenic (As) Cadmium (Cd) Mercury (Hg) Lead (Pb) Expiration Date Period 5 Years

Specification: PRD.0.ZQ5.10000066566

Specification

Colorless to Yellow Liquid 1.442 - 1.446 Conforms to Structure > 97.5 % < 500 APHA < 3 ppm < 1 ppm < 1 ppm < 10 ppm

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.