

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Monoclonal Anti-polyHistidine-Alkaline Phosphatase antibody produced in mouse - clone HIS-1, purified immunoglobulin, buffered aqueous glycerol solution

Product Number:

**A5588**

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Almost Colorless to Light Yellow
IgG Concentration (OD 280nm)	0.8 - 2.5 mg/ml
WB-Cell Line/Tissue Extract	E.Coli Exp PolyHistidine Tagged Protein
WB Titer	1:2000-1:4000
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000015098

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.