

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lipodisq® - Styrene:Maleic Anhydride Copolymer 3:1

Product Number:

L9045

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
IR Spectrum	Conforms to Structure
Specific Gravity	Conforms
1.05-1.34 at 25 Deg C	
Registered Trademark	Conforms
Lipodisq is a registered trademark of Malvern Cosmeceutics Limited	
Composition	98 - 100 %
by size exclusion chromatography(SEC)	

Specification: PRD.1.ZQ5.10000037311

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

