

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platelet-Derived Growth Factor-BB human – PDGF-BB, recombinant, expressed in E. coli, lyophilized powder, suitable for cell culture

Product Number:

P3201

CAS Number:

105186-99-0

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Form)	Lyophilized Powder
Bioactivity	Conforms
EC50 cell proliferation with Balb 3T3 cells	
1.0-3.0 ng/ml	
Purity	Conforms
Purity by SDS-PAGE and HPLC, > / = 98%	
Endotoxin Assay	< or equal to 1 Eu per ug

Specification: PRD.1.ZQ5.10000093443

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.