

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver bis(trifluoromethanesulfonyl)imide - 97%

Product Number:

668001

CAS Number:


189114-61-2

Formula:

C₂AgF₆NO₄S₂

Formula Weight:

388.01 g/mol

**TEST****Specification**

Appearance (Color)

White to Beige to Orange

Appearance (Form)

Powder or Crystals or Chunk(s) or

Powder with Chunk(s) or Crystals with Chunk(s)

Infrared spectrum

Conforms to Structure

ICP Major Analysis

Confirmed

Confirms Ag and S Components

Titration by KSCN

26.8 - 28.8 %

% Ag

Specification: PRD.1.ZQ5.10000012623

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.