

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Nitrate Reductase (NAD[P]H) from *Aspergillus niger* - lyophilized powder, ≥ 300 units/g solid

Product Number:

N7265

CAS Number:

9029-27-0

MDL:

MFCD00131725

Storage Temperature:

-20 °C

TEST**Specification**

units/G solid

 ≥ 300

Enzymatic Activity: One unit will reduce 1.0 micromole of nitrate per minute in the presence of beta-NADPH at pH 7.5, 25 deg C.

Expiration Dated Product

Specification: PRD.3.ZQ5.10000008700

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.