

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Betamethasone - ≥98%

Product Number: B7005 CAS Number: 378-44-9

 MDL:
 MFCD00062969

 Formula:
 C22H29FO5

 Formula Weight:
 392.46 g/mol

HO HOCH

TEST Specification

Appearance (Color) White

Appearance (Form) Powder EMM 15.0 - 15.6

at Lambda Max 238 to 240 nm, in MeOH

Wavelength 238 - 240 nm

Solubility (Color) Colorless
Solubility (Turbidity) Clear

20 mg + 1 mL, EtOH:CHCl3 (1:1)

Purity (HPLC) \geq 98 %

Specification: PRD.0.ZQ5.10000045645

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.