


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,2-Bis[(2S,5S)-2,5-diisopropylphospholano]benzene(1,5-cyclooctadiene)rhodium(I) tetrafluoroborate

Product Number: **698334**
CAS Number: 1000393-36-1
MDL: MFCD07369040
Formula: C₃₄H₅₆BF₄P₂Rh
Formula Weight: 716.47 g/mol



TEST	Specification
Appearance (Color) Yellow to Red to Brown	Conforms to Requirements
Appearance (Form) Powder, Crystals or Crystalline Powder	Conforms to Requirements
Proton NMR spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Carbon	54.7 - 59.3 %
Optical Rotation c = 0.8%; Chloroform	78.0 - 86.0 deg

Specification: PRD.0.ZQ5.10000052769

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.