

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Diethylene glycol monoethyl ether acetate - 99%

Product Number:

537527

CAS Number:

112-15-2

MDL:


MFCD00041928

Formula:

C₈H₁₆O₄

Formula Weight:

176.21 g/mol

**TEST****Specification**

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 0.0500 %
Specific Gravity at 20 Degrees Celsius	1.008 - 1.013
Acidity (as Acetic Acid)	≤ 0.0200 %
Color Test	≤ 10 APHA
Assay	≥ 99.00 %

Specification: PRD.0.ZQ5.10000017743

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.