

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Proteose Peptone – Enzymatic hydrolysate

Product Number: **P0431**

TEST	Specification
------	---------------

Appearance (Color)	Faint Yellow to Yellow and Faint Beige to Beige and Faint Brown to Light Brown
Appearance (Form)	Powder
Solubility (Color)	Faint Yellow to Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
at 1% w/v in H ₂ O	
Nitrogen	10.0 - 15.5 %
Loss on Drying	≤ 4 %
pH	6.6 - 7.6
(1% w/v in H ₂ O)	

Specification: PRD.1.ZQ5.10000023243

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.