

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Purity (HPLC)

Ammonium

Salicylic acid - BioXtra, ≥99.0%

Product Number: S5922 CAS Number: 69-72-7

MDL: MFCD00002439

Formula: C7H6O3 Formula Weight: 138.12 g/mol

TEST Specification

Appearance (Color) White to Off-White Appearance (Form) Powder or Crystals Infrared spectrum Conforms to Structure Solubility (Turbidity) Clear 1M in Ethanol at 20 Degrees Celsius Solubility (Color) Colorless Residue on ignition (Ash) < 0.01 % Insoluble matter 0.1 % Chloride Content 0.001 % Sulfate 0.003 % < 0.005 % Aluminum (AI) Calcium (Ca) 0.001 % Copper (Cu) 0.0005 % Iron (Fe) 0.0005 % < 0.005 % Potassium (K) Magnesium (Mg) 0.0005 % Sodium (Na) 0.01 % Phosphorus (P) 0.015 % < 0.001 % Lead (Pb) Zinc (Zn) 0.0005 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

99.0 %

< 0.05 %



3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Salicylic acid - BioXtra, ≥99.0%

Product Number: \$5922 CAS Number: 69-72-7

MDL: MFCD00002439

Formula: C7H6O3

Formula Weight: 138.12 g/mol

TEST Specification

Recommended Retest Period

4 Years

Specification: PRD.1.ZQ5.10000058791