

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lithium hydroxide - reagent grade, 98%

Product Number:**442410**

LiOH

CAS Number:

1310-65-2

MDL:

MFCD00011095

Formula:

HLiO

Formula Weight:

23.95 g/mol

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Conforms to Requirements

Flakes, Chips, Granules, Pellets, or Chunks

Titration by HCL

97.5 - 102.5 %

ICP Major Analysis

Confirmed

Confirms Lithium Component

Specification: PRD.1.ZQ5.10000004921

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.