


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Bromoquinoxaline - 97%Product Number: **752339**
CAS Number: 36856-91-4Formula: C₈H₅BrN₂
Formula Weight: 209.04 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Light Grey to Very Dark Grey and Light Brown to Dark Brown
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	44.3 - 47.6 %
Nitrogen	12.9 - 13.9 %

Specification: PRD.0.ZQ5.10000022497

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.