

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

NaK Silica gel (K₂Na) Stage I - extent of labeling: ~35 wt. % Na-K loading

Product Number:

660140

CAS Number:

11135-81-2

K₂Na

TEST	Specification
------	---------------

Appearance (Color)

Black

Appearance (Form)

Powder

ICP Major Analysis

Confirmed

Confirms Sodium, Potassium and Silicon Components

Active Material Loading

30 - 35 % wt

K₂Na (wt%)

Specification: PRD.0.ZQ5.10000081968

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.