

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Sea salts

Product Number: S9883

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
38 mg/ml, Water	
Chloride (Cl)	Conforms
19000-20000 mg/L	
Sodium	Conforms
10700-11000mg/L	
Sulfate	Conforms
2660 mg/L	
Potassium	Conforms
300-400 mg/L	
Calcium	Conforms
400 mg/L	
Carbonate	Conforms
140-200mg/L	
Boron	Conforms
5.6 mg/L	
Magnesium (Mg)	Conforms
1320 mg/L	
Strontium (Sr)	Conforms
8.8 mg/L	
Insoluble Matter	≤ 0.05 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Sea saltsProduct Number: **S9883**

TEST	Specification
------	---------------

Note	Conforms
Typical Ion Concentration of sea water. Also contains < 0.5 mg/L other trace elements.	

Specification: PRD.1.ZQ5.10000031869

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.