


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Hydroxyanthranilic acid - 97%

Product Number: 148776
CAS Number: 548-93-6
MDL: MFCD00007700
Formula: C₇H₇NO₃
Formula Weight: 153.14 g/mol



TEST

Specification

Appearance (Color)	Off White to Light Yellow to Light Orange
Appearance (Form)	Powder or Crystals
Solubility (Color)	Yellow to Brown-Yellow and Light Orange to Orange
Solubility (Turbidity)	Clear
10 mg/ml, 1N HCl	
Infrared spectrum	Conforms to Structure
Carbon	53.0 - 56.8 %
Nitrogen	8.8 - 9.5 %
Purity (HPLC)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000034480

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.