

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethylamine solution - 2.0 M in methanol

Product Number: **395064**



MDL: MFCD00008160

Formula: C₂H₇N

Formula Weight: 45.08 g/mol

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Titration with NaOH

1.90 - 2.20 M

Specification: PRD.1.ZQ5.10000065814

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.