


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Isoamyltriphenylphosphonium bromide - 98%

Product Number: 347779
CAS Number: 28322-40-9
MDL: MFCD00031594
Formula: C₂₃H₂₆P · Br
Formula Weight: 413.34 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Light Orange
Appearance (Form)	Crystals or Powder or Granules or Chunks or Crystals with Chunks
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 0.1 %
Titration by AgNO ₃	≥ 97.5 %

Specification: PRD.1.ZQ5.10000034985

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

