

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Carbon nanotube array, multi-walled – vertically aligned on Si substrate, H 0.5 mm

Product Number:

901051

CAS Number:

308068-56-6



Formula:

C

TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Solid
ICP Major Analysis	Confirmed
Confirms Silicon Component	
Length	Confirmed
(Height: 0.5mm +/- 0.1mm)	
Diameter	Confirmed
10-40 nm	
Substrate	Confirmed
Dimension 20x20 mm	
Carbon Content	≥ 98 %

Specification: PRD.2.ZQ5.10000086712

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.