

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Copper(II) sulfate pentahydrate - 99.999% trace metals basis

Product Number:**203165**CuSO₄ • 5H₂O

CAS Number:

7758-99-8

MDL:

MFCD00149681

Formula:

CuO4S · 5H2O

Formula Weight:

249.69 g/mol

TEST	Specification
Appearance (Color)	Blue
Appearance (Form) Crystals and/or Chunks	Conforms to Requirements
ICP Major Analysis Confirms Copper and Sulfur Components	Confirmed
Purity 99.999% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.0.ZQ5.10000085620

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.