


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Benzoic anhydride - ≥95%

Product Number: 385980
CAS Number: 93-97-0
MDL: MFCD00003073
Formula: C₁₄H₁₀O₃
Formula Weight: 226.23 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystals or Granules or Chunks
Infrared spectrum	Conforms to Structure
Morpholine Titration by HCl	≥ 95.0 %

Specification: PRD.1.ZQ5.10000020864

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.