


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Iron(II) sulfide - technical grade

Product Number: 268704
CAS Number: 1317-37-9
MDL: MFCD00011013
Formula: FeS
Formula Weight: 87.91 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Grey To Black	
Appearance (Form)	Powder or Chunks
ICP Major Analysis	Confirmed
Confirms Iron and Sulfur Components	
Titration with Na ₂ S ₂ O ₃	60.0 - 72.0 %
% Iron	

Specification: PRD.1.ZQ5.10000033282

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.