

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Gelatin solution – Type B, 2% in H₂O, tissue culture grade, BioReagent, suitable for cell culture**Product Number:** G1393

MDL: MFCD00081638

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Turbidity) at Room Temperature	Slightly Hazy
Appearance (Color) at Room Temperature	Faint Yellow
Appearance (Form) at Room Temperature	Solution
Protein by Biuret	14 - 22 mg/ml
Bioburden Report results	
Endotoxin Level Report result	
Cell Culture Test Cell attachment: Greater than or equal to 80% of control	≥ 80 %
Cell Culture Test Cell Toxicity: Greater than or equal to 80% of control	≥ 80 %

Specification: PRD.6.ZQ5.10000012613

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.