

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Albumin human – recombinant, expressed in rice, lyophilized powder, suitable for cell culture, low endotoxin,
≥96% (SDS-PAGE)

Product Number:

A9731

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

White to Light Yellow to Light Brown

Appearance (Form)

Powder

Solubility (Color)

Very Faint Yellow to Yellow to Brown-Yellow

Solubility (Turbidity)

Clear to Hazy

50 mg/mL, H₂O

Water (by Karl Fischer)

≤ 12 %

Cell Culture Test

Pass

Endotoxin Level

≤ 1 EU/mg

Note

Cellastim is a trademark of InVitria

Purity by SDS-PAGE

≥ 96 %

Specification: PRD.5.ZQ5.10000020100

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

