

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
rac-BINAP Pd G4

Product Number: 900340

Formula: C₅₈H₄₇NO₃P₂PdS
Formula Weight: 1,006.45 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow to Brown
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
³¹ P NMR Spectrum	Conforms to Structure
Carbon Content	62.0 - 76.5 %
Nitrogen Content	
Specification: Report Only	
Sulfur Content	
Specification: Report Only	
ICP Major Analysis	Confirmed
Confirms Palladium Component	

Specification: PRD.0.ZQ5.10000079682

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.