

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum – wire, diam. 1.0 mm, 99.9% trace metals basis

Product Number:**327492**

CAS Number:

7440-06-4

MDL:

MFCD00011179

Formula:

Pt

Formula Weight:

195.08 g/mol

Pt

TEST	Specification
Appearance (Color) Silver	Conforms
Appearance (Form) Wire	Conforms
Diameter 1.0 mm	Conforms
ICP Major Analysis Confirms Platinum Component	Confirmed
Purity 99.9% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 1500.0 ppm

Specification: PRD.1.ZQ5.10000003619

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.