

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(ethylene glycol) dimethacrylate – average M_n 6,000, contains 1000 ppm 4-methoxyphenol as inhibitor

Product Number:


687537

CAS Number:

25852-47-5

MDL:

MFCD00081877



Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Inhibitor	Confirmed
1000 ppm 4-methoxyphenol	
Molecular Number	5400 - 6600
Poly Dispersity MW/MN	

Specification: PRD.0.ZQ5.10000093548

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.