

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phytase from wheat – ≥ 0.01 unit/mg solid

Product Number:

P1259

CAS Number:

9001-89-2

MDL:

MFCD00131910

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/mg solid	≥ 0.01
----------------	-------------

Enzymatic Activity

Unit Definition: One unit will liberate 1 micromole of

inorganic phosphorus from 1.5 MM phytate per minute at pH 5.15 at 55 deg C

Specification: PRD.0.ZQ5.10000081692

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.