

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3,4-Ethylenedioxy pyrrole - 2 % (w/v) in THF

Product Number:

648310

CAS Number:

169616-17-5

MDL:

MFCD07369727

Formula:


C₆H₇NO₂

Formula Weight:

125.13 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

Dark Brown to Very Dark Brown and Black

Appearance (Form)

Liquid

Proton NMR spectra

Conforms to Structure

Specification: PRD.0.ZQ5.10000087252

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.