

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(tetrahydrofuran) - average $M_n \sim 1,000$

Product Number: **345296**
 CAS Number: 25190-06-1
 MDL: MFCD00148879
 Formula: $H(OCH_2CH_2CH_2CH_2)_nOH$



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Colorless	
Appearance (Form)	Conforms to Requirements
Waxy Solid to Viscous Liquid	
Infrared spectrum	Conforms to Structure
Molecular Number	950 - 1050
Hydroxyl Value (mg KOH/g)	107.0 - 118.0
Color Test	≤ 40 APHA

Specification: PRD.1.ZQ5.10000003707

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.