

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cobalt(II) 2-ethylhexanoate solution - 65 wt. % in mineral spirits

Product Number: **444545**

MDL: MFCD00072632
Formula: C₁₆H₃₀CoO₄
Formula Weight: 345.34 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Dark Blue to Very Dark Blue and Dark Purple to Very Dark Purple and Blue-Purple
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Complexometric EDTA % Co	11.0 - 13.0 %

Specification: PRD.0.ZQ5.10000038710

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.