

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Fluorophosphoric acid solution - 70 wt. % in H₂O

Product Number:


324744

MDL:

MFCD00011334

Formula Weight:

99.99 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Brown

Appearance (Form)

Liquid

Titration with NaOH

60.0 - 80.0

Specification: PRD.1.ZQ5.1000005101

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.