


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Ethylphenol - 99%

Product Number: **E44205**
CAS Number: 123-07-9
MDL: MFCD00002393
Formula: C₈H₁₀O
Formula Weight: 122.16 g/mol

**TEST****Specification**

Appearance (Color)
Appearance (Form)
Crystals or Solid
Infrared spectrum
GC (area %)

White to Yellow and White to Orange
Conforms to Requirements
Conforms to Structure
≥ 98.5 %

Specification: PRD.2.ZQ5.10000005600

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.