

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Rhenium powder, -100 mesh, ≥99.9% trace metals basis

Product Number:**267279**

CAS Number:

7440-15-5

MDL:

MFCD00011195

Formula:

Re

Formula Weight:

186.21 g/mol

Re

TEST**Specification**

Appearance (Color)

Grey

Conforms to Requirements

Appearance (Form)

Powder

X-Ray Diffraction

Conforms to Structure

Trace Metal Analysis

≤1,000.0 ppm

Purity

≥99.9% Based on Trace Metals Analysis

Conforms

Size

-100 Mesh

Conforms

Remarks:

Specification Date : 11/17/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.