

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:

Ammonium molybdate - 99.98% trace metals basis

Product Number:**277908**

CAS Number:

13106-76-8



MDL:

MFCD00011596

Formula:

H8MoN2O4

Formula Weight:

196.01 g/mol

TEST**Specification**

Appearance (Color)

White

Appearance (Form)

Powder or Crystals or Chunk(s) or Solid
or Crystalline Chunk(s) or Powder with Chunk(s) or Crystals with
Chunk(s)

Gravimetric Analysis

47.2 - 50.7 %

% Mo With Ammonium Molybdate

ICP Major Analysis

Confirmed

Confirms Molybdenum Component

Purity

Meets Requirements

99.98% Based On Trace Metals Analysis

Trace Metal Analysis

≤ 250.0 ppm

Specification: PRD.1.ZQ5.10000039191

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.