

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium pentafluoroethyltrifluoroborate - 95%

Product Number: **747157**
CAS Number: 476639-90-4

Formula: C₂BF₈K
Formula Weight: 225.92 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	White to Faint Yellow and Faint Beige and Faint Brown
Appearance (Form)	Powder or Crystals
Fluorine NMR Spectra	Conforms to Structure
Boron NMR Spectrum	Conforms to Structure
Carbon	10.0 - 11.2 %
ICP Major Analysis	Confirmed
Confirms K Component	

Specification: PRD.1.ZQ5.10000041075

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.