

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lead(II) selenide - 99.99%

Product Number: 254231
CAS Number: 12069-00-0
MDL: MFCD00016273
Formula: PbSe
Formula Weight: 286.16 g/mol

PbSe

TEST**Specification**

Appearance (Color)	Conforms to Requirements
Appearance (Form)	Powder or Granules or Chunk(s)
ICP: Confirms Pb & Se Componen	Confirmed
Purity	Conforms
99.99% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 200.0 ppm

Specification: PRD.0.ZQ5.10000015847

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.