


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Methylresorcinol - 97%

Product Number: 657581
CAS Number: 496-73-1
MDL: MFCD00016452
Formula: C₇H₈O₂
Formula Weight: 124.14 g/mol



TEST

Specification

Appearance (Color)

White to Yellow and Faint Beige to Light

Appearance (Form)

Beige and Faint Red

Solid, Fibers, Powder, Crystals, Crystalline Powder and/or Chunks

Conforms to Requirements

Infrared spectrum

Conforms to Structure

GC (area %)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000026213

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.